

Project Title: Pediatric Dose Calculator

SAT 4650 FINAL PROJECT

BY IFUNANYA EZEUMEH

Project description

- Plan: To develop a Python program that will help calculate compounding doses for pediatric populations.
- Modules: Functions, Database & GUI.
- Dataflow/ major functionalities:
 - Users will use entry widget to input patient relevant information.
 - Formulas to be used to calculate the doses will be defined with functions and stored to a database.
 - The calculated doses will be outputted in an info dialog box.
- Input: Relevant patient information, and medication name
- Output: Calculated compounding dose

Major Application of The Pediatric Dosage Calculator

Children of the same age grow at different rates. For children in the same age group, there are wide variations in their weight . Consequently, dosages are usually calculated using their weight. It is important to ensure that every child is getting the optimum and safe dose required for them to receive quality healthcare.

Long Term Goals

To expand the functions of the calculator to be able to calculate doses for other special populations. Thereafter, create an Application Programming Interface (API) that can be embedded into a Pharmacy Health Information Technology (HIT) solution which will help Pharmacists quickly determine doses required for many special populations.

```

1 import tkinter
2 from tkinter import *
3 from tkinter.ttk import *
4 import mysql.connector
5 from tkinter import messagebox as tkMessageBox
6 from DosageCalculator_db_fxns import *
7
8 conn = mysql.connector.connect(
9     host= "localhost",
10    user= "root",
11    passwd= "password")
12
13 db_cursor = conn.cursor()
14
15 class MyGUI:
16
17     def __init__(self):
18         # create the main window
19         self.main_window = tkinter.Tk()
20         self.main_window.title("Pediatric Dosage Calculator")
21         self.main_window.configure(background='aliceblue')
22         self.calculated_dose_value = StringVar()
23
24         # create a menu
25         menu = Menu(self.main_window)
26         self.main_window.config(menu=menu)
27         filemenu = Menu(menu)
28         menu.add_cascade(label="File", menu=filemenu)
29         filemenu.add_command(label="New Patient", command=callback)
30         filemenu.add_command(label="Open Recent", command=callback)
31         filemenu.add_separator()
32         filemenu.add_command(label="Exit Calculator", command=callback)
33         helpmenu = Menu(menu)
34         menu.add_cascade(label="Help", menu=helpmenu)
35         helpmenu.add_command(label="About Pediatric Dosage Calculator", command=callback)
36
37         # create frames to group widgets
38         self.top_frame = tkinter.Frame(self.main_window, bg= "aliceblue")
39         self.mid_frame = tkinter.Frame(self.main_window, bg= "aliceblue")
40         self.bottom_frame = tkinter.Frame(self.main_window, bg= "aliceblue")
41
42         # Top frame
43         # create widgets
44         self.note_label = tkinter.Label(self.top_frame, text="Note: If child is weighed in pounds (lbs) "

```

```

49 self.weight_entry = tkinter.Entry(self.top_frame, font="bold", width = 30, justify="left")
50 self.weight_entry.grid(row=1, column=1)
51 self.weight_label = tkinter.Label(self.top_frame, bg= "aliceblue", font= ("bold", 14), text="Enter Weight (kg): ", justify="left")
52 self.weight_label.grid(row=1, column=0)
53
54 self.dosing_rule_entry = tkinter.Entry(self.top_frame, font="bold", width = 30, justify="left")
55 self.dosing_rule_entry.grid(row=2, column=1)
56 self.dosing_rule_label = tkinter.Label(self.top_frame, bg= "aliceblue", font= ("bold", 14), text="Enter Dosing Rule (mg/kg/day): ", justify="left")
57 self.dosing_rule_label.grid(row=2, column=0)
58
59 self.frequency_entry = tkinter.Entry(self.top_frame, font= "bold", width = 30, justify="left")
60 self.frequency_entry.grid(row=3, column=1)
61 self.frequency_label = tkinter.Label(self.top_frame, bg= "aliceblue", font= ("bold", 14), text="Enter Frequency (day): ", justify="left")
62 self.frequency_label.grid(row=3, column=0)
63
64 self.strength_entry = tkinter.Entry(self.top_frame, font= "bold", width = 30, justify="left")
65 self.strength_entry.grid(row=4, column=1)
66 self.strength_label = tkinter.Label(self.top_frame, text="Enter Dosage strength in mg/5ml: ", bg= "aliceblue", font= ("bold", 14), justify="left")
67 self.strength_label.grid(row=4, column=0)
68
69 self.disease_label = tkinter.Label(self.top_frame, text= "Indication List: ", bg= "aliceblue", width = 30, font= ("bold", 14), justify="left")
70 self.disease_label.grid(row=5, column=0)
71 data1 = ("Acute Otitis Media", "Allergy Relief", "Pain Relief")
72 self.disease_combobox = Combobox(self.top_frame, values = data1, width= 28, font= "bold")
73 self.disease_combobox.grid(row=5, column=1)
74
75 self.medication_label = tkinter.Label(self.top_frame, bg= "aliceblue", text="Medication List: ", width=30, font= ("bold", 14), justify="left")
76 self.medication_label.grid(row=6, column=0)
77 data2 = ("Amoxicillin", "Diphenylhramine", "Ibuprofen", "Tynelol")
78 self.medication_combobox = Combobox(self.top_frame, values = data2, font= "bold", width=28 )
79 self.medication_combobox.grid(row=6, column=1)
80
81 # Middle frame
82 # create widgets for mid frame
83 self.cal_button = tkinter.Button(self.mid_frame, bg= "alice blue", text = "Calculate", font= ("bold", 14), command=self.calculateDose)
84 self.cal_button.grid(row=9, column=0, columnspan=1, pady=10, padx=10, ipadx=66)
85
86 # create store button for mid frame
87 self.save_button = tkinter.Button(self.mid_frame, bg= "alice blue", font= ("bold", 14), text = "Save", command= self.save)
88 self.save_button.grid(row=9, column=1, columnspan=1, pady=10, padx=10, ipadx=50)
89
90 # create reset button for mid frame
91 self.reset_button = tkinter.Button(self.mid_frame, bg= "alice blue", font= ("bold", 14), text = "Reset", command= self.reset)
92 self.reset_button.grid(row=9, column=2, columnspan=1, pady=10, padx=10, ipadx=66)

```



```

103         .....# call the label widget's pack method
104         .....self.top_frame.pack()
105         .....self.mid_frame.pack()
106         .....self.bottom_frame.pack()
107
108         .....tkinter.mainloop()
109
110         .....# To reset the calculator ...
111         .....def reset(self):
112         .....    conn = mysql.connector.connect(
113         .....        host= "localhost",
114         .....        user= "root",
115         .....        passwd= "password",
116         .....        database = "dosageCalculator_db")
117         .....    db_cursor = conn.cursor()
118         .....
119         .....# To clear texboxes
120         .....    self.weight_entry.delete (0, END)
121         .....    self.strength_entry.delete (0, END)
122         .....    self.frequency_entry.delete (0, END)
123         .....    self.dosing_rule_entry.delete (0, END)
124         .....    self.disease_combobox.set("")
125         .....    self.medication_combobox.set("") .....
126         .....    conn.commit()
127         .....    conn.close()
128
129         .....#To store to DB and clear fields.
130         .....def save(self):
131         .....    if self.weight_entry.get() == "" or self.strength_entry.get() == "" or self.frequency_entry.get() == "" or self.dosing_rule_entry.get() == "":
132         .....        .....tkMessageBox.showerror('Error!', 'Enter correct details')
133         .....    else:
134         .....        .....conn = mysql.connector.connect(
135         .....            host= "localhost",
136         .....            user= "root",
137         .....            passwd= "password",
138         .....            database = "dosageCalculator_db")
139         .....        .....db_cursor = conn.cursor()
140
141         .....# To insert data into a particular column
142         .....        sql_wt = "Insert into Weight_Info (weight_kg) values (%s)"
143         .....        db_cursor.execute(sql_wt, (self.weight_entry.get(),))
144         .....        conn.commit()
145         .....        sql_str = "Insert into Strength_Info (dose_strength_mg_per_5ml) values (%s)"

```

```

216         frequency = db_cursor.execute(sql_fq)
217         result_frequency = db_cursor.fetchall()
218         for frequency in result_frequency:
219             input_freq = (float(frequency[0]))
220         else:
221             sql_fq = "Insert into Frequency_Info (frequency_per_day) values (%s)"
222
223             frequency = db_cursor.execute(sql_fq, (self.frequency_entry.get(),))
224             input_freq = float(self.frequency_entry.get())
225             conn.commit()
226             print(input_freq)
227             print(type(input_freq))
228             conn.commit()
229
230         dose_rule = float(self.dosing_rule_entry.get())
231         if dose_rule == "":
232             sql_dr = "select dosing_rule_mg_per_kg_per_day from Dosing_Rule_Info where class='A'"
233             dose_rule = db_cursor.execute(sql_dr)
234             result_dose_rule = db_cursor.fetchall()
235             for dose_rule in result_dose_rule:
236                 input_dose = (float(dose_rule[0]))
237             else:
238                 sql_dr = "Insert into Dosing_Rule_Info (dosing_rule_mg_per_kg_per_day) values (%s)"
239
240                 dose_rule = db_cursor.execute(sql_fq, (self.dosing_rule_entry.get(),))
241                 input_dose = float(self.dosing_rule_entry.get())
242                 conn.commit()
243                 print(input_dose)
244                 print(type(input_dose))
245                 conn.commit()
246
247         mg_per_day = (input_weight * input_dose)
248         mg_per_dose = mg_per_day / input_freq
249         dose_to_ml = mg_per_dose / input_strength
250         dose_per_ml = round(dose_to_ml, 2)
251
252         print(dose_per_ml)
253         self.calculated_dose_value.set("Recommended dose is: " + str(dose_per_ml) + " mL to be taken " + str(input_freq) + " time(s) daily.")
254         conn.commit()
255         conn.close()
256
257 my_gui = MyGUI()

```


Pediatric Dosage Calculator


Note: If child is weighed in pounds (lbs), divide by 2.2 to obtain weight in kilograms(kg). For example, 22 lbs / 2.2 = 10kg.

Enter Weight (kg):

Enter Dosing Rule (mg/kg/day):

Enter Frequency (day):

Enter Dosage strength in mg/5ml:

Indication List: 

Medication List: 